

Mitochondrial respiration in permeabilized fibres: needle biopsies from horse skeletal muscle



Oxygraph-2k Workshop Protocol (IOC44),
December 2007, Schroecken, Austria.

Lemieux H¹, Votion DM²,
Doerrier C³, Gnaiger E^{1,3}

¹D. Swarovski Research Lab, Dept Visceral, Transplant and Thoracic Surgery,
Medical Univ Innsbruck, Austria, www.mitofit.org

²Dept of Clinical Sciences, Fac of Veterinary Medicine and Equine
European Centre of Mont-le-Soie, Univ of Liège, Sart Tilman, 4000 Liège, Belgium

³Oroboros Instruments, Schöpfstr 18, A-6020 Innsbruck, Austria,
Email: instruments@orooboros.at; www.orooboros.at

Section

1. Introduction	1
2. The SUIT protocol and respiratory states.....	2
3. References	4

1. Introduction

Methodological and conceptual features of High-Resolution Respirometry (HRR) are illustrated in an experiment with permeabilized fibres in the Oroboros Oxygraph-2k (O2k). A mitochondrial substrate-uncoupler-inhibitor titration (SUIT) protocol with manual titrations was applied to study physiologically relevant maximum mitochondrial respiratory capacity and coupling/pathway control. Experiments were carried out

by participants of an O2k-Workshop in December 2007 (IOC44; Schroecken, Austria; Votion et al 2012).

2. The SUIT protocol and respiratory states




Figure 1. Oxygen concentration ($[\mu\text{M}]$ blue line) and oxygen flux per mg wet weight of muscle ($[\text{pmol}\cdot\text{s}^{-1}\cdot\text{mg}^{-1}]$ red line) in permeabilized fibres from horse skeletal muscle.




[\3.O2k-Procedures\MiPNet12.23_FibreRespiration\MiPNet12.23_Pfi_2007-03-21 P1-01.DLD](#)

2.1. The O2k-Demo experiment

Permeabilized fibres from horse skeletal muscle (*Triceps brachii*) were prepared (Pesta and Gnaiger 2012) and incubated at 37 °C in the Oxygraph-2k, with 2 mL of mitochondrial respiration medium (MiR05 or MiR06 [MiPNet14.13]).

2.2. SUIT events, marks, and respiratory states



In the SUIT protocol (Fig. 1) a sequence of respiratory states is induced experimentally by stepwise titrations (Events, E). As a consequence of the titrated compounds, respiration reaches a new steady-state, and a mark (M) is set for numerical evaluation of the corresponding respiratory state.

1GM;2D;2c;3S;4U;5Rot;6Ama

E 1G,1M

10 mM glutamate & 2 mM malate was added to the chambers before adding the fibres (1.5 to 2.5 mg wet weight), resting state.

M 1GM

GM_L or GM_L: LEAK state with type N substrates, **N_L:** NADH-linked substrates glutamate&malate (type N; CI-linked pathway to Q). Non-phosphorylating resting state (LEAK state L; in the absence of ADP; no adenylates).

- E **2D** After titration of 2.5 mM ADP (D), flux increases to active respiration.
- M **2D** **GM_P or GM_P: OXPHOS capacity with type N substrates, N_P:** Respiratory capacity in the active coupled OXPHOS state (with saturating [ADP], *P*).
- E **2c** 10 μM cytochrome *c* is added as a test for the intactness of the mitochondrial outer membrane (mtOM).
- M **2c** **GMc_P or GMc_P: Cytochrome *c* test for quality control:** Addition of cytochrome *c* yields a test for integrity of the outer mitochondrial membrane (loss of cytochrome *c* would be indicated by a stimulation of respiration).
- E **3S** Respiration is further stimulated by succinate (10 mM).
- M **3S** **GMS_P or GMS_P: OXPHOS capacity with type NS substrates (CI&II-linked pathway to Q), NS_P:** Respiratory stimulation by further addition of succinate (*S*) to type N substrates, with convergent electron flow through CI and CII at the Q-junction, in the coupled state, as an estimate of OXPHOS capacity with reconstitution of the TCA cycle ([Gnaiger 2009](#)); *P*, OXPHOS capacity of the NS-pathway.
- E **4U*** Subsequently, FCCP (uncoupler, U) is titrated in steps of 0.125 μM, to test for a possible increase of noncoupled flux compared to state GMS_P.
- M **4U*** **GMS_E or GMS_E: Electron transfer (ET) capacity with type NS substrates, NS_E:** Uncoupling by FCCP titration (avoiding inhibition by high FCCP concentrations), as a test for limitation of OXPHOS-capacity by the phosphorylation system (ANT, ATP synthase, phosphate transporter) relative to ET-capacity (ADP activated, noncoupled state, *E*, ETS capacity or the NS-pathway).
- E **5Rot** Inhibition of CI by rotenone (0.5 μM).
- M **5Rot** **S_E or S_E: ET capacity with type S (CII) substrate:** S-pathway ET-capacity after blocking CI with rotenone.
- E **6Ama** Inhibition of CIII by Antimycin A (Ama; 2.5 μM) or myxothiazole.
- M **6Ama** **Residual oxygen consumption (ROX)** due to oxidative side reactions, estimated after addition of Antimycin A (inhibitor of CIII). Respiration of uncoupled mitochondria might be inhibited slightly further by cyanide (KCN; 1 μM). ROX is subtracted from oxygen flux as a baseline for all respiratory states, to obtain mitochondrial respiration.




Figure 2 A: Mitochondrial O₂ flux corrected for ROX. **B:** Flux control ratios normalized to ETS capacity.



Excel demo file:

[\3.O2k-Procedures\MiPNet12.23_FibreRespiration\MiPNet12.23_Pfi_O2k-Analysis.xls](#)

3. References

- Gnaiger E (2009) Capacity of oxidative phosphorylation in human skeletal muscle. New perspectives of mitochondrial physiology. *Int J Biochem Cell Biol* 41:1837–45. [»Bioblast link«](#)
- Pesta D, Gnaiger E (2012) High-resolution respirometry. OXPHOS protocols for human cells and permeabilized fibres from small biopsies of human muscle. *Methods Mol Biol* 810:25-58. [»Bioblast link«](#)
- Votion DM, Gnaiger E, Lemieux H, Mouithys-Mickalad A, Serteyn D (2012) Physical fitness and mitochondrial respiratory capacity in horse skeletal muscle. *PLoS ONE* 7(4):e34890. [»Bioblast link«](#)